


Comfortable to breath
Hypoallergenic
Latex free
Soft elastic ear loops
Plastic wire nose support
Outer layer colored non woven | black
Middle layer meltblown BFE≥98%
Inner layer non woven



Surgical BFE ≥98% TYPE II | EN 14683:2019



20 pcs

How to wear surgical mask

- 1 Clean your hands with soap and hand sanitizer
- 2 Hold the mask by the ear loops and place a loop around each ear
- 3 Mold or pinch the stiff edge to the shape of your nose
- 4 Pull the bottom of the mask over your mouth and chin
- 5 Avoid touching the front of the mask when wearing


How to remove surgical mask

- 1 Clean your hands with soap and water or hand sanitizer
- 2 Avoid touching the front of the mask. Only touch the ear loops
- 3 Hold both of the ear loops and gently lift and remove the mask
- 4 Throw the mask in the trash
- 5 Clean your hands with soap and water or hand sanitizer

Hellenic Medical Supplies GP 6th km Xanthi-Lefkopetra, 67100 Xanthi, Greece, Europe T +30 25410 32300 E info@masks.gr www.masks.gr



STORE UNDER COOL, CLEAN AND DRY CONDITIONS



Made in Clean Room, Class D

Surgical BFE ≥98% TYPE II | EN 14683:2019

Surgical face mask is intended to be worn by medical personnel during surgical or other medical procedures to limit the transmission of infective agents from staff to patients during surgical procedures and other medical settings with similar requirements.

20 pcs

- | | | | | |
|---------------------------|----------------------------|---------------------------|------------------------|------------------------|
| GR/ES Χειρουργικές μάσκες | EE Kirurgiline näomask | IE Mascanna máinliachta | PL Maski chirurgiczne | SK Chirurgicalké masky |
| BG Хирургически маски | ES Mascarillas quirúrgicas | IT Mascherine chirurgiche | PT Máscaras cirúrgicas | TR Cerrahi maske |
| CRO Kirurške maske | FI Kirurgiset naamiot | LT Chiruginės kaukės | RO Măști chirurgicale | AZ Cerrahi maskalar |
| CZ Chirurgicalké roušky | FR Masques chirurgicaux | LV Maskas kirurģiskās | RU Хирургические маски | |
| DK Operations masker | GB/US Surgical masks | MT Maski kirurģiċi | SE Kirurgiska masker | |
| DE Chirurgische masken | HU Sebészeti maszkok | NL Chirurgische maskers | SI Kirurške maske | |

RECYCLABLE PAPER PACKAGE